

Multifunctional Skin Care Platform 8 in 1

GUARDIAN Of Your Youth !!

- ▶ Deep Cleaning
- ▶ Skin Tightening
- ▶ Acne Treatment
- ▶ Moisturizing
- ▶ Skin Repairing
- ▶ Essence Import
- ▶ Skin Analysis






1. Dull complexion




2. Large pores



3. Wrinkles



4. Dry Skin



5. Acne



6. Sagging skin



ONE MACHINE

Solves All Skin Problems


HI-SKIN multifunctional skin care instrument is a comprehensive eight kinds of latest technologies, which can solve almost all facial skin problems. It is the best choice for skin clinics and beauty salons. The instrument is equipped with many of the most commonly used treatment functions, including cleaning, hydrating, lifting and firming, acne, sedation, dead skin removal, skin analysis, etc.






HYDROGEN OXYGEN HYDRO DERMABRASION HANDLE


Open pores, soften the stratum corneum, deeply cleanse the skin, deeply moisturize, remove black-heads, improve the dark yellow skin, and make the skin clean, white and supple.



Solution is sprayed through the inlet of nozzle tip. It rotates clockwise.



The waste and remaining solution are sucked to container through inlets



Solution are supplied to skin surface through nozzle tip

A Variety Of Treatment Heads Are Equipped, Users Can Choose The Clean Treatment Head They Need According To The Skin Condition



Hard turbine hydro cleaning head



Microdermabrasion head



Soft silicone hydro cleaning head



Pipe cleaning head

COOL RF AND COLD HAMMER HANDLE

Cool RF handle uses an intelligent epidermal cooling system. By pre-cooling the epidermis, most of the RF energy is pushed to the deep layers of the skin. While protecting the epidermis from burns, the heat energy is concentrated in the deep layers of the skin to stimulate collagen regeneration. The two functions of RF and Cold hammer can be used separately. Turn off the RF energy, The handle can also be used as an Cold hammer alone to shrink pores, calm the epidermis, and eliminate local redness.



Polymer Atomizing Spray Handle

Combined with liquid products, the liquid is converted into mist-like small molecules by high-pressure atomization, and then quickly poured into the dermis by high-pressure spraying to achieve a deep hydrating effect







PLASMA BIOPOLAR RF HANDLE

Plasma technology can be effectively used for anti-inflammatory and sterilization, remove mites, and treat acne. Bipolar radio frequency is used to improve fine lines at the corners of the eyes, remove bags and dark circles under the eyes.

NANO MESO MICROWAVE INTRODUCING HANDLE

Adopt nano-microwave technology to open the skin channel with micro-wave stimulation, effectively enhance the skin's absorption of nutrients, and solve the skin absorption problem painlessly and noninvasively





HD SKIN DETECTION PROBE

The high-definition skin detection probe uses a high-definition camera with 2,000,000 pixels and a resolution of 1280 * 1024, which can clearly observe the true condition of the skin and every pore.



ULTRASONIC VIBRATION SKIN CLEANING HANDLE

Remove dead skin by high-frequency vibration waves, remove black-heads, soften skin keratinocytes, and deeply cleanse skin dirt.



SCOPE OF TREATMENT

Deep Cleaning

Remove Blackheads, Mites
Clear, Remove grease dust,
Remove Dead Skin

Skin Tightening

Face Lifting, Eye Lifting, Wrinkle
removal

Skin Rejuvenation

Skin Moisturizing, Shrink Pores,
Essence Import, Skin Whitening,
Acne Treatment

Skin Analysis

Analyze skin conditions and
compare skin conditions before
and after treatment

BEFORE AND AFTER



SPECIFICATION

Power	600w
RF frequency	3MHz
Screen	10.4inch touch screen
Pixel	20million
Probe of dermabrasion	10 pics
Vacuum	0-90Kpa
Flow Capacity	30 L Max
Voltage	100V-240V 50Hz-60Hz